

F-DECALIN PLUS



More Information:
www.fluoron.de

Biocompatible

Physiologically and chemically inert

Adjusted density of 1.81 g/cm³

F-DECALIN^{PLUS} is a mixture of Perfluorodecalin and F6H8[®]. Its reduced density promotes "gentle" reattachment of the retina.

- Intraoperative adjuvant for retinal reattachment.
- The addition of F6H8[®] allows remaining droplets to build a homogenous tamponade medium with the ensuing silicone oil tamponade.
- The reduced density allows tractions to remain visible to the surgeon.

G-80205 F-DECALIN^{PLUS}
Vial 5 ml, sterile

G-80207 F-DECALIN^{PLUS}
Vial 7 ml, sterile