

INNOVATION
MADE BY GEUDER®

NEW



PURE EFFICIENCY PHACO

Greatest Efficiency for
Smallest Incisions

AUFFARTH
MICS DOUBLE CROSS ACTION
CAPSULORHEXIS FORCEPS



G-32926 AUFFARTH
MICS DOUBLE CROSS ACTION
CAPSULORHEXIS FORCEPS
flat handle, cross-action,
for incision size 1.6 mm,
total length 100 mm



Prof. Gerd U. Auffarth, MD,
Clinical Director
University Eye Hospital Heidelberg

More freedom of movement through the smallest clear cornea incision (down to 1.6 mm):

- The narrow design guarantees great freedom of movement in the anterior chamber.
- The hand-crafted and precise forceps tips facilitate safe handling of the capsulorhexis.
- The overlying design of the forceps tips prevents the iris from getting caught in them, especially in cases of intraoperative floppy iris syndrome (IFIS).

Distributed by:

IQ Medical Pty Ltd | 2/86 Mary Street, Unley SA 5061

Phone (08) 8357 8022 | Email sales@iqmedical.com.au

GEUDER AG
Hertzstrasse 4
69126 Heidelberg
Germany

Phone: +49 6221 3066
Fax: +49 6221 303122
info@geuder.de
www.geuder.de

Geuder®
Precision made in Germany