


SURGILŪM

See what you've been missing

ROBO-MARKER™

The world's first self-leveling corneal marker with pre-inked, sterile, disposable tips and an integrated fixation light.






Used for Toric IOL Placement

Ink marks using the RoboMarker's dry-ink formula last up to two hours and maintain clear, crisp lines. No autoclave wait and single step marking means a more efficient surgical day.

ROBO-MARKER

Profiled as one of the top eight innovations at ASCRS 2014, Surgilūm's state of the art RoboMarker redefines corneal marking for astigmatic surgery. The RoboMarker places corneal marks within 1 degree of accuracy that last up to two hours. By redesigning the process of marking, the single-use RoboTip offers surgeons a safe and efficient method for achieving long-lasting precise marks, an efficient surgical day, and brilliant patient outcomes.



Used for Femtosecond Laser Alignment

Premium results depend upon accurate marks. The RoboMarker's self-leveling high precision ball bearing based pendular mechanism ensures that your Femtosecond Laser is aligned properly for optimal surgical results.



ROBOTIP

The disposable companion to RoboMarker.

Sterile

No more waiting for the autoclave

Pre-inked

**Our unique gentian violet
ink lasts up to 2 hours**

Disposable

**The single-use RoboTip makes
marking your patients as
easy as putting on a new tip**

You are a Premium Surgeon.

You invest exorbitantly to gather premium pre-op data.

Your patient trusts you for a premium result. And now...

Your personal premium RoboMarker. Let it take you there.

Efficient • Accurate • Elegant

ROBO  **MARKER**™



VIEW THE **ROBOMARKER** IN ACTION



ALAN W. BROWN, MD

Inventor of the RoboMarker™

Founder, Surgilūm

MEDICAL ADVISORY BOARD



KARL STONECIPHER, MD

“The greatest thing I like about the RoboMarker™ is its ease of use and the fact that the marks last for over an hour. It really has improved my outcomes and my patient flow.”



MITCHELL A. JACKSON, MD

“RoboMarker™ is the most accurate way to mark for intended toric axis alignment prior to toric IOL and/or astigmatic corneal incision surgery due to its automatic cyclotorsion adjustment.”



GEORGE O. WARING IV, MD

“The RoboMarker™ is an advancement in the surgical planning of astigmatism correction. The intelligent balance takes the guesswork out of your alignment and is intuitive to use.”



WILLIAM TRATTLER, MD

“The RoboMarker™ is very easy to use in the Preop area, and once marked, it makes aligning toric IOLs very easy.”



R. DAX HAWKINS, MD

“When marking the axis for surgery, I need something that’s accurate, quick and easy. The RoboMarker™ is all of these in a package that fits in the chest pocket of my scrubs.”



RON KRUEGER, MD

“When astigmatism correction is needed to make a difference in outcomes, the Robomarker™ makes a difference in the axis alignment.”

See what you've been missing



Distributed by:

IQ Medical Pty Ltd

2/86 Mary Street, Unley SA 5061

Phone (08) 8357 8022

Email sales@iqmedical.com.au

Web www.iqmedical.com.au