

F-OCTANE PLUS



More Information:
www.fluoron.de

Biocompatible

Physiologically and chemically inert

Reduced density of 1.72 g/cm³

F-OCTANE PLUS is a mixture of Perfluorooctane and F6H8®. Its reduced density promotes "gentle" reattachment of the retina.

- Intraoperative adjuvant for retinal reattachment.
- The reduced density allows tractions to remain visible to the surgeon.
- The four-times greater vapor pressure allows remaining droplets to vaporize faster than perfluorodecalin droplets, providing quicker relief for the retina.

G-80405 F-OCTANE PLUS
Vial 5 ml, sterile

G-80407 F-OCTANE PLUS
Vial 7 ml, sterile