

INNOVATION  
MADE BY GEUDER®

NEW

# PURE EFFICIENCY PHACO




Greatest Efficiency for  
Smallest Incisions

AUFFARTH  
MICS DOUBLE CROSS ACTION  
CAPSULORHEXIS FORCEPS



G-32926 AUFFARTH  
MICS DOUBLE CROSS ACTION  
CAPSULORHEXIS FORCEPS

flat handle, cross-action,  
for incision size 1.6 mm,  
total length 100 mm



Prof. Gerd U. Auffarth, MD,  
Clinical Director  
University Eye Hospital Heidelberg

### More freedom of movement through the smallest clear cornea incision (down to 1.6 mm):

- The narrow design guarantees great freedom of movement in the anterior chamber.
- The hand-crafted and precise forceps tips facilitate safe handling of the capsulorhexis.
- The overlying design of the forceps tips prevents the iris from getting caught in them, especially in cases of intraoperative floppy iris syndrome (IFIS).

**Geuder**®  
Precision made in Germany

GEUDER AG  
Hertzstrasse 4  
69126 Heidelberg  
Germany

Phone: +49 6221 3066  
Fax.: +49 6221 303122  
info@geuder.de  
www.geuder.de