


## One Piece Equi-aspheric Foldable Acrylic

Optic Size : 5.75mm  
 Optic Type : Equiconic  
 Length : 12.00mm  
 Haptic Style: Modified C  
 Angulation : 0 Degrees  
 Construction : 1 Piece  
 Optic Material : 26% Hydrophilic Acrylic  
 A Constant : 118.0  
 A/C Depth : 5.10mm



- Precision manufactured to greater than ISO standard with mid range tolerance of 0.125D
- Square edge technology

**Power Range :** -5.0 to +36.0  
**Standard Range :** +10.0D to +30.0D  
**Special Order :** -5.0 to +9.0 and +31.0 to +36.0  
**0.25 Diopter Steps** from +18.0 to +25.0  
**0.5 Diopter Steps** from +10.5 to +18.0  
 and +25.5 to +30.0

Please note : If you are implanting Softec HD greater than +30.0D, a 2.7 cartridge system must be used.

Modified A Constant and A/C Depth values for IOL Master users shown below are strictly guidelines for the calculation of implant power. IQ Medical suggests surgeons develop their own values based on technique, measuring equipment and desired post-op results.

### IOL Master

SRKII/SRKT	A Constant	118.2		
Holl 2	ACD	4.907		
Holl 1	SF	1.276		
Hoff Q	ACD	5.083		
Contact	A Constant	118.0		
Haigis	a0 = 0.92	a1 = 0.40	a2 = 0.10	
Immersion	A Constant	same as IOL Master		

When using IOL Master please ensure .25 box is selected in lens set up. The appropriate surgical techniques are the responsibility of the respective surgeon. They must assess the appropriateness of the relevant procedure based on their education and experience.

**For further information and ordering Softec HD & HDY please contact IQ Medical**

on Phone: (08) 8357 8022, FAX: (08) 8357 8011

or Email: [sales@iqmedical.com.au](mailto:sales@iqmedical.com.au)

Please allow 3-5 working days for delivery.