

Manufacturers of specialist surgical instruments



SINGLE USE PRODUCT CATALOGUE

Manufacturers of specialist surgical instruments

We are **Surgitrac® Instruments** a multinational manufacturer and distributor of specialist surgical instruments for use in Ophthalmology.

Our manufacturing site is ISO 13485:2016 certified.

Surgitrac Instruments UK is ISO 13485: 2016 and EN ISO 13485:2016 certified.

We offer:

- Quality assured CE marked products, manufactured to internationally recognised standards.
- Transparency in our supply chain, no forced, or child labour is used during our manufacturing.
- A continuous supply of products, due to our capacity to manufacture millions of instruments each year.

We are specialists in single-use instruments

Ophthalmic instruments arrive sterile, and ready to use.

We offer any combination of instruments assembled as sterile, sets and procedure packs.

Procedure packs include a selection of instruments and associated items such as drapes, swabs and gallipots.

Instruments sets include two or more items, sterile and in a convenient wrap or tray ready to use.

Packs that we offer include:

Cataract | DMEK basic | DMEK complete | Eye bank | Cross Linking | Lasek | Lasik | Intravitreal injection (IVT)
Vitreoretinal (VR) | Lacrimal wash-out | Chalazion | Sub-tenon anaesthesia.

Sets that we offer include:

Cataract | Glaucoma | Pterygium | Biopsy & minor procedure | Eye lid surgery | Plastic surgery
Sub tenon anaesthesia | Strabismus | Retina of Prematurity (ROP) | Eye bank

Non sterile bulk wholesale

Our instruments are available for bulk wholesale. Instruments are supplied packed and ready for sterilisation. We can offer private label manufacturing through our OEM business. Precision instruments featuring your brand can be supplied as single-use sterile or as non-sterile units.

Contact us info@surgitrac.com more information.

Using this catalogue

This catalogue is designed to provide specific information about our products and their technical specifications.

Product details

Product code number – Single use S codes.



■ Single Use Instrument Code

Capsulorhexis forceps cross action - SC33

Length: 11.5 cm
Handle length: 82 mm
50° angles platform



Description details product specifications such as handle length.

Product images are not shown at actual size.

The technical illustration specifies details such as tip measurements.

Detailed technical files for instruments are available at request.

CONTENTS:

Page 2: About

Page 4: User guide

Page 5: Contents

Page 6: Forceps

Page 18: Manipulators and Choppers

Page 24: Speculums

Page 30: Scissors

Page 36: General ophthalmic

Page 43: Corneal and Refractive

Page 52: Measurers and Markers

Page 56: Curettes and Clamps

Page 59: Needle holders

Page 61: Lacrimal

Page 64: Index

Forceps


Capsulorhexis | Fixation | General | Lens | Tying



Single Use Instrument Code


Capsulorhexis forceps cross action - SC33

Length: 11.5 cm
 Handle length: 82 mm
 50° angled platform tips 0.3 mm




Capsulorhexis forceps cross action no tips - SC154

Length: 11.7cm
 Handle length: 85 mm




Capsulorhexis forceps curved - SC32

Length: 10.5 mm
 Handle length: 75 mm
 45° angled platform with tips 0.3 mm




Capsulorhexis forceps round bodied - SC34

Length: 10.5 cm
 Handle length: 77 mm
 45° angled platform with tips 0.3 mm



Capsulorhexis forceps straight - SC30

Length: 80 mm
 Handle length: 56 mm
 50° angled platform tips 0.5 mm




Capsulorhexis forceps straight • SC31

Length: 10.5 cm
 Handle length: 75 mm
 45° angled
 12 mm platform tips 0.3 mm



CAPSULORHEXIS

**MICS Capsulorhexis forceps cross action** • SC33M

Length: 11.5 cm
 Handle length: 84 mm
 50° angled
 2.2 mm incision



CAPSULORHEXIS

**Capsulorhexis forceps curved** • SC29

Length: 10.5 cm
 Handle length: 76 mm
 45° angled platform tips 0.3 mm 2.4 mm
 incision



CAPSULORHEXIS

**MICS Capsulorhexis forceps curved** • SC29M

Length: 10.5 cm
 Handle length: 77 cm
 45° angle
 12 mm platform
 0.3 mm tips
 2.2 mm incision



CAPSULORHEXIS




MICS Capsulorhexis forceps round bodied curved • SC34M

Length: 10.5 cm
 Handle length: 77 mm
 50° angle 0.2 mm tip
 12 mm platform
 2.2 mm incision



CAPSULORHEXIS



MICS Capsulorhexis forceps straight • SCM01

Length: 10.5 cm
 Handle length: 76 mm
 45° angled 0.3 mm tips, with 11 mm platform 2.2 mm incision



CAPSULORHEXIS



MICS Capsulorhexis forceps straight tipped • SC133M

Length: 10.5 cm
 Handle length: 76 mm



CAPSULORHEXIS



MICS Capsulorhexis forceps round bodied straight tipped • SC134M

Length: 10.6 cm
 Handle length: 76 mm



CAPSULORHEXIS




Adsons dissection forceps - SC59

Length: 12.5 cm
 Handle length: 97 mm
 1 into 2 teeth
 28 mm platform




FIXATION

**Bishops forceps - SC56**

Length: 10.5 cm
 Handle length: 75 mm
 1 into 2 teeth
 5 mm platform




FIXATION

**Bonn forceps - SC38**

Length: 85 mm
 Handle length: 57 mm
 1 into 2 teeth 5 mm platform
 Short handle
 5 mm platform




FIXATION

**Bonn forceps - SC129**

Length: 10.5 mm
 Handle length: 76 mm
 1 into 2 teeth 5 mm platform



FIXATION

**Colibri toothed forceps - SC22**

Length: 11 cm
 Handle length: 72 mm
 55° angled
 1 into 2 teeth
 4 mm platform



FIXATION




Colibri tying forceps - SC20

Length: 85 mm
Handle length: 59 mm
55° angle forked tip
Short handle
4 mm platform



FIXATION




Fine claw forceps - SC115

Length: 10 cm
Handle length: 68 mm
2 into 1 teeth
32 mm platform to tip



FIXATION



Hoskins forceps - SC55

Length: 85 mm
Handle length: 56 mm
Notched forceps
5 mm platform



FIXATION



Iris forceps - SC116

Length: 10.5 cm
Handle length: 67 mm
1 into 2 teeth 5 mm platform



FIXATION




Iris forceps pointed - SC128

Length: 11 cm
Handle length: 69 mm
Serrated 15 mm platform



FIXATION




Jayles forceps • SC117

Length: 10.5 cm
 Handle length: 67 mm
 1 into 2 teeth 5 mm platform




FIXATION

**Straight notched forceps** • SC49


Length: 11 cm
 Handle length: 76 mm
 Notched forceps
 5 mm platform




FIXATION

**St. Martins forceps** • SC54

Length: 84 mm
 Handle length: 56 mm
 1 into 2 teeth
 6 mm platform




FIXATION

**Straight toothed forceps fine** • SC27

Length: 10.5 cm
 Handle length: 78 mm
 1 into 2 teeth 5 mm platform




FIXATION




Adsons forceps • SC119

Length: 12.5 mm
 Handle length: 97 mm
 14 mm serrated platform




Adsons forceps fine serrated • SC113

Length: 12.5 cm
 Handle length: 88 mm
 serrated 14 mm tying platform




Denis sclerectomy forceps • SC104

Length: 10.5 cm
 Handle length: 75 mm
 30 mm platforms



Epilation - Cilia forceps • SC58

Length: 73 mm
 Platform 50 mm




Epilation forceps - SC127

Length: 88 mm
18 mm platform




GENERAL

**Jewellers forceps** - SC60

Length: 11.5 cm
Handle length: 82 mm
0.1 mm fine pointed ends
32 mm platform




GENERAL

**Moorfields forceps** - SA12

Length: 11 cm
Handle length: 68 mm
Serrated 15 mm platform




GENERAL

**Moorfields fine forceps** - SA17

Length: 11 cm
Handle length: 71 mm
Serrated 15 mm platform




GENERAL




Mosquito forceps - SG02

Length: 13 cm
22 mm serrated clamp area




IOL removal forceps - SC112

Length: 82 mm
Handle length: 54 mm
35° angled serrated tips, with 13 mm platform




Lens folding forceps - SC51

Length: 10.5 cm
Handle length: 78 mm
45° angled, 1 cm platform



Lens implantation forceps - SC52

Length: 10.5 cm
Handle length: 76 mm
40° angled platform
10 mm platform




Curved notched forceps • SC50

Length: 10 cm
 Handle length: 76 mm
 Curved 5 mm tying platforms




TYING

**Kelman Mc Pherson forceps** • SC73

Length: 83 mm
 Handle length: 57 mm
 50° angled 7 mm tying platforms




TYING

**Mc Pherson tying forceps** • SC72

Length: 10.5 cm
 Handle length: 76 mm
 45° angled 5 mm tying platforms




TYING

**Straight tying forceps** • SC19

Length: 85 mm
 Handle length: 56 mm
 5 mm tying platforms




TYING




Straight tying forceps - SC71

5 mm platforms
Length: 10.5 cm
Handle length: 75 mm
5 mm tying platform




Suture removal forceps curved - SC120

Length: 10.5 cm
Handle length: 76 mm
Pointed tips, 5 mm tying platform



Troutman tying forceps curved - SC105

Length: 10.5 cm
Handle length: 76 mm
5 mm tying platform




Manipulators and Choppers



Single Use Instrument Code


Chopper and dialler - SC39

Length: 12 cm
 Handle length: 76 mm
 45° angled ends, one rounded tip and one 1 mm Y shaped end
 Round bodied
 Double ended




Chopper and kuglen horizontal double ended - SC95

Length: 12.7 cm
 Handle length: 70 mm
 40° angled paddle shaped tip
 Round bodied




Denis cleaver - SC97

Length: 11 cm
 Handle length: 90 mm



Drysdale micro manipulator - SC89

Length: 12 cm
 Handle length: 90 mm
 Round bodied



Drysdale nucleus spatula - SC47

Length: 11 cm
 Handle length: 85 mm
 45° angled paddle shaped tip
 Round bodied




Govan manipulator - SC100


60° angle pointed tip
Round bodied

**Govan manipulator and Phaco chopper - SC102**

Length: 12 cm
Handle length: 78 mm
60° pointed end and 90° blunt end
Round bodied
Double ended

**Greens chopper - SC45**

Length: 11.5 cm
Handle length: 95 mm
30° angled top curved rounded hook
Round bodied

**Iris repositor - SC63**

Length: 13 cm
Handle length: 98 mm


**Jones chopper - SC90**

Length: 12 cm
Handle length: 95 mm
45° angled pointed tip
Round bodied




Koch chopper - SC44

Length: 12 cm
 Handle length: 95 mm
 45° angled 1 mm rounded tip
 Round bodied



Koch phaco spatula - SC143

Length: 12 cm
 Handle length: 96 mm
 Slight curved spatula
 Round bodied



Kuglen hook - SG15

Length: 11.5 cm
 Handle length: 95 mm
 45° angled blunt, flat hooked tip
 Round bodied



Le Buisson manipulator - SC106

Length: 11.5 cm
 Handle length: 90 mm
 Round bodied



Lens dialer - SC46

Length: 12 cm
 Handle length: 95 mm
 35° angled
 Round bodied




Lester hook • SC146


Length: 11.5 cm
 Handle length: 10 cm
 Round bodied

**Mushroom nucleus manipulator** • SC48


Length: 12 cm
 Handle length: 95 mm
 45° angled mushroom shaped tip
 Round bodied

**Phaco chopper and Mushroom manipulator double ended** • SC147

Length: 12 cm
 Handle length: 70 mm
 45° chopper
 Round bodied
 Double ended

**Phaco chopper** • SC43

Length: 12 cm
 Handle length: 95 mm
 40° angled round end
 Round bodied



Phaco chopper claw - SC42

Length: 12 cm
Handle length: 95 mm
40° angled tip, hooked tip
Round bodied




Pisacano rotator and Neuham manipulator double end - SC101

Length: 10 cm
Handle length: 75 mm
60° spatula end and 30° round end
Round bodied
Double ended



Sinskey hook and Govan manipulator double ended - SC144

Length: 12 cm
Handle length: 76 mm
60° manipulator and 45° hooked end
Round bodied
Double ended



Y shaped Nucleus rotator - SC41

Length: 12 cm
Handle length: 95 mm
40° angled tip Y shaped end
Round bodied



Speculums

Wire | Mechanical | Paediatric




■ Single Use Instrument Code

Broad gauge wire speculum - SC01

Length: 50 mm
Width: 40 mm




WIRE




Heavy gauge wire speculum open loop - SC124

Length: 72 mm
Overall width: 58 mm



WIRE



Kratz Barraquer closed loop - SC09

Length: 40 mm
Width: 42 mm




WIRE




Kratz Barraquer closed loop speculum - SC132

Length: 43 mm
50° angle at base




WIRE



Kratz Barraquer open loop wire speculum - SC11


Length: 47 mm
 Width: 31 mm
 50° angle at base




WIRE

**Kratz open loop wire speculum (blue) - SC107**


Length: 44 mm
 Overall width: 30 mm
 50° angle at base



WIRE

**Kratz solid bladed wire speculum - SC05**


WIRE

**Kratz closed loop speculum - SC13**

Length: 32 mm
 50° angle at base




WIRE

**Open loop wire speculum - SC150**

Length: 38 mm
 Width: 38 mm



WIRE




Speculum open loop • SC155

Length: 45 mm
Overall width: 30 mm




WIRE



Temporal speculum • SC12

Length: 50 mm
Width: 42 mm



WIRE




Wide angled closed loop speculum • SC07

Length: 44 mm
Width: 42 mm



WIRE




Bangerter adult speculum • SC03

Length: 50 mm
sliding screw adjustable
Blade length: 15 mm



WIRE




Bladed adjustable speculum - SR10

Length: 75 mm
Blade length: 15 mm




MECHANICAL

**Clarke speculum - SC126**

Length: 80 mm
Open loop: 12 mm




MECHANICAL

**Kratz screw adjustable wire speculum open loop - SR11**

Length: 78 mm
Blade length 15 mm




MECHANICAL

**Lieberman speculum stainless steel - SC121**

Length: 80 mm
Blade length 14 mm




MECHANICAL

**Solid bladed screw adjustable speculum - SC08**

Length: 85 mm
Blade length: 14 mm



MECHANICAL




Kratz paediatric closed loop 10 mm speculum - SC10

Length: 29 mm
Width: 30 mm
10 mm



PAEDIATRIC



Neonate paediatric speculum - SC103

Length: 25 mm
Width: 20 mm
4 mm closed loop
80° angle at base



PAEDIATRIC



Neonate 4 mm speculum - SC15

Length: 25 mm
Width: 31 mm
20-30 mm adjustable
3 mm blade



PAEDIATRIC




Paediatric 6 mm speculum - SC14

Length: 25 mm
Width: 35 mm
6 mm blade



PAEDIATRIC



Scissors

Spring Handled | General




Single Use Instrument Code

Castroviejo corneal scissor right and left • SC82

Length: 10 cm
 Handle length: 75 mm
 ■ SC83 Castroviejo corneal scissor left
 8 mm blade



SPRING HANDLED




Castroviejo scissors curved and blunt • SC110

Length: 11 cm
 Handle: 65 cm
 11 mm blade



SPRING HANDLED




Castroviejo scissors straight • SC111

Length: 10 cm
 Handle length: 63 mm
 Straight tipped
 60 mm blade



SPRING HANDLED




Ong iris scissor • SC80

Length: 88 mm
 Handle: 52 mm
 10 mm blade



SPRING HANDLED




Osher scissors • SC109

Length: 95 mm
 Handle length: 61 mm
 Fine serrated cutting platform
 9 mm blade



SPRING HANDLED




Stallard corneal scissors • SA13

Length: 12 cm
 Handle length: 67 mm
 Fine pointed ends
 40 mm blade




SPRING HANDLED

**Vannas scissors** • SC81

Length: 83 mm
 Handle length: 52 mm
 8 mm blade




SPRING HANDLED

**Vannas scissors curved** • SC108

Length: 10 cm
 Handle length: 62 mm
 Sharp pointed end
 9 mm blade



SPRING HANDLED

**Westcott conjunctival scissors** • SA11

Length: 11.5 cm
 Handle length: 79 cm
 12 mm blade



SPRING HANDLED




Wescott scissors dolphin nose • SA14

Length: 11.5 cm
Handle length: 67 mm
10 mm blade



SPRING HANDLED




Westcott dolphin nose scissors blunt tipped • SA15

Length: 11.6 cm
Handle length: 72 mm



SPRING HANDLED



Dressing scissors • SG35

Length: 11.5 cm
Handle length: 26 mm
30 mm blade



GENERAL



Enucleation scissors • SG34

Length: 25 mm
6 mm blade



GENERAL




General long handle scissors blunt end • SG14

Length: 13 cm
Blunt round end
40 mm blade




GENERAL

**General scissors blunt end • SG06**

Length: 89 mm
Blunt end
24 mm blade




GENERAL

**General scissors sharp pointed • SG01**

Length: 88 mm
Sharp pointed end
22 mm blade




GENERAL

**Iris scissor • SC96**

Length: 11 cm
22 mm blade



GENERAL




Metzenbaum scissors curved pointed ▪ SG30

Length: 11 cm
35 mm blade



GENERAL



Strabismus scissors ▪ SG07

Length: 11.5 cm
Blunt end
38 mm blade



GENERAL




Stevens scissors curved ▪ SG23

Length: 11.5 cm
38 mm blade



GENERAL




Stevens tenotomy scissors pointed ▪ SG31

Length: 11.6 cm
Tip width: 1.9mm



GENERAL




General ophthalmic



■ Single Use Instrument Code


Ball tip muscle hook • SG32

Length: 10.5 cm
Handle length: 76 mm




Bulldog clip • SG92

Length: 40 mm
Handle length: 25 mm
16 mm clipping platform not on website




Chavasse hook • SC91

Length: 13 cm
Handle length: 98 mm
90° angle, curved hooked tip




Desmarres retractor 18mm wide • SG20

Length: 14 cm
Handle length: 100 mm




Eyelid plate • SG17


Length: 90 mm
27 mm end
17 mm end

**Gillies hook small** • SG28

Length: 12 cm
Handle length: 90 mm

**Guthrie hook** • SG37


Length: 12 cm
Handle length: 90 mm
Two pronged hook

**Halstead mosquito forceps straight** • SG33

Length: 12.5 cm


**Jameson strabismus hook 5 mm** • SG10

Length: 11 cm
Handle length: 83 mm
60° hook with rounded bulbous tip
• SG11 extended version 8 mm hook



Mosquito forceps curved • SG22

Length: 12.5 cm
22 mm serrated clamp area



Muscle hook with hole • SG18

Length: 13.8 cm
Handle length: 95 mm



Scleral depressor bulbous end • SC23

Length: 13 cm
Handle length: 92 mm
2.5 mm 10° angled bulbous end




Scleral depressor double ended • SC25

Length: 13.5 cm
Handle length: 60 mm
6 mm, 2 mm hammer end,
2.5 mm 10° angled bulbous end
Double ended




Scleral depressor hammer shape • SC24

Length: 12.5 cm
Handle length: 92 mm
6 mm hammer, 2 mm ends




Thornton fixation ring small 13 mm • SG12


Length: 13 cm
 Handle length: 91 mm
 swivel head 11 mm cut out


**Thornton fixation ring standard 16 mm • SG13**

Length: 13.5 cm
 Handle length: 91 mm
 swivel head 14 mm cut out


**Trabeculectomy punch • SG03**

Length: 12 cm
 Handle length: 90 mm





Corneal and Refractive




Single Use Instrument Code

Alcohol well various sizes - SR41

Length: 12.5 cm / Handle length: 92 mm
Liquid holding well for alcohol application / Round bodied


- SR41 9 mm
- SR42 8 mm
- SR43 7 mm

**Bernheim hook DSAEK - SR109**

Length: 12 cm
Handle length: 95 mm

**Bladed adjustable speculum - SR10**

Length: 46 mm
Blade length: 15 mm
Blade radius: 2.25mm

**Corneal epithelium spatula round - SR100**

Length: 11 cm
Handle length: 85 mm




Corneal marker 8 line • SR72


Length: 11 cm
Handle length: 82 mm

**CK Lasik marker** • SR40

Length: 10 cm
Handle length: 76 mm
facilitate marking of the cornea
Round bodied

**Corneal ring pusher** • SR04

Length 11 cm
Handle length 91 mm

**Corneal ring pocketing hook** • SR05

Length 11 cm
Handle length 91 mm




Corneal ring pocket spreader - SR06

Length 10 cm
Handle length 90 mm




Corneal ring forceps Keraring - SR07

Length 10cm
Handle length 77 mm



Corneal ring forceps INTACS - SR08

Length 8 cm
Handle length 53 mm




Epi peeler rectangle shaped - SR61

Length 11 cm
Handle length 95 mm

**Epi peeler D shaped** - SR60


Length 11 cm
Handle length 95 mm

**Epi peeler** - SR51


Length: 11 cm
Handle length: 95 mm
Manipulate epithelium during surgery

**Implantation forceps** - SR91

Length: 117 mm
Handle length: 70 mm


**Intra Lasik flap separator** - SR34

Length: 12 cm
Round bodied




Kocaba DMEK forceps - SR102

Length: 86 mm
Handle length: 56 mm




Kratz screw adjustable open loop speculum - SR11

Length: 78 mm
Blade length: 15 mm




Lasik corneal marker - SR21

Length: 12 cm
Handle length: 91 mm
Y shaped forked marker



Lasik femto separator and elevator double ended - SR02

Length: 115 mm
Handle length: 70 mm




Lasik flap elevator - SR31


Length: 12 cm
 Handle length: 95 mm
 45° angled, 12 mm tip with 0.1 mm rounded end

**Lasik flap forceps - SR14**


Length: 10 cm
 Handle length: 76 mm

**Lasik flap forceps - SR71**

Length: 10 cm
 Handle length: 77 mm
 40° angled platform, circular 1 mm tips


**Lasik flap separator and elevator double end - SR103**

Length: 12 cm
 Handle length: 70 mm
 45° and 650° ends
 Double ended



Lenticule extraction forceps • SR01

Length: 115 mm
 Handle length: 81 mm




Paton spatula and spoon double ended • SR108

Length: 151 mm
 Handle length: 71 mm




Reinstein Lenticule separator elevator double ended • SR104

Length 11 cm
 Double ended
 Handle length 70 mm




Reinstein Lenticule separator • SR70

Length: 110 mm
 Handle length: 90 mm




S Marker - SR101


Length: 11 cm
 Handle length: 90 mm
 45° angle 2 mm marking end


**Sinskey and reverse Sinskey hook** - SR03

Length 12 cm
 Handle length 70 mm

**Trephine various sizes** - SR105

- SR110 16 mm
- SR105 17 mm
- SR106 18 mm
- SR107 19 mm






Measurers and Markers



Single Use Instrument Code


Braunstein fixed caliper 1.8 mm thickness - SG04

Length: 68 mm
 3.5 mm 4 mm measure
 (Available in other thicknesses of metal, please contact us for details)



Castroviejo adjustable caliper - SG05

Length: 85 mm
 20 mm measure




Corneal marker capsulorhexis guide - SC130

Length: 12 cm
 Handle length: 92 mm
 6 mm well
 Round bodied




Heleveston ruler - SG93


Length: 11.7 cm
Handle length: 89 mm


**Mendez Gauge** - SC65

Length: 12 cm
Handle length: 80 mm
180° circular marker

**Scleral marker round bodied double ended** - SG19

4.5 and 3.5 mm marking ends
Length: 11 cm
Handle: 60 mm
Round bodied
Double ended





Curettes and Clamps




■ Single Use Instrument Code

Curette various size

Length: 11 cm
 Handle length: 80 mm
 ■ CR02 Curette medium, 2.5 mm diameter chalazion curette
 ■ CR03 Curette large, 4 mm oval shaped chalazion curette



CURETTES




Curette double ended - CR05

Length: 13cm
 Handle length: 65 mm
 chalazion curette 2.5 mm diameter end and oval end



CLAMPS




Chalazion clamp small - CS01

Length: 11.5 cm
 Handle length: 82 mm
 12 mm diameter plate



CLAMPS




Chalazion clamp medium - CS03

Length: 90 mm
 Handle length: 54 mm
 19 mm diameter plate



CLAMPS




Chalazion clamp large - CS02

Length: 85 mm
 Handle length: 53 mm
 21 mm diameter oval plate



CLAMPS




Needle holders



Single Use Instrument Code


Barraquer needle holder curved • SC142

Length: 11.5 cm
 Handle length: 82 mm
 12 mm diameter plate
 Round bodied




Castroviejo needle holder locking • SC118

Length: 14 cm




Needle holder blunt lock • SC62

Length: 13.5 cm
 Handle length: 82 mm




Needle holder curved • SC114

Length: 10 cm
 Handle length: 75 mm



Suture needle holder curved • SC61

Length: 13.5 mm
 Handle length: 86 mm
 6 mm jaw



Lacrimal



■ Single Use Instrument Code

Bowmans Probe size 00-0 double ended • PD05


Length: 50 mm
Handle length: 30 mm

**Bowmans Probe size 000-00 double ended** • PD06

Length: 50 mm
Handle length: 30 mm

**Lacrimal probe** • PD08

Length: 95 mm
Handle length: 15 mm
0.5 mm pointed end


**Punctal dilator blunt** • PD02

Length: 95 mm
Handle length: 80 mm
0.4 mm pointed sharp end
Round bodied




Punctal dilator double ended - PD03


Length: 11.5 cm
 Handle length: 85 mm
 0.4 and 0.3 mm sharp pointed ends
 Round bodied

**Punctal dilator sharp - PD01**

Length: 95 mm
 Handle length: 78 mm
 0.2 mm pointed sharp end
 Round bodied


**Punctal dilator short - PD04**

Length: 50 mm
 Handle length: 30 mm
 0.3 mm sharp pointed ends
 Round bodied

**Punctal dilator & lacrimal probe double ended - PD07**

Length: 14 cm
 Handle length: 70 mm
 curved and pointed end





CAPSULORHEXIS CAPSULORHEXIS KAPSULORHEXIS - Pg7-9

| | |
|---------------|--|
| SC33 | CAPSULORHEXIS FORCEPS CROSS ACTION |
| SC154 | CAPSULORHEXIS FORCEPS CROSS ACTION |
| SC32 | CAPSULORHEXIS FORCEPS CURVED |
| SC34 | CAPSULORHEXIS FORCEPS ROUND BODIED |
| SC30 | CAPSULORHEXIS FORCEPS STRAIGHT |
| SC3 | CAPSULORHEXIS FORCEPS STRAIGHT |
| SC33M | MICS CAPSULORHEXIS FORCEPS CROSS ACTION |
| SC29 | CAPSULORHEXIS FORCEPS CURVED |
| SC29M | MICS CAPSULORHEXIS FORCEPS CURVED |
| SC34M | MICS CAPSULORHEXIS FORCEPS ROUND BODIED CURVED |
| SCM01 | MICS FORCEPS STRAIGHT |
| SC133M | MICS CAPSULORHEXIS FORCEPS STRAIGHT TIPPED |
| SC134M | MICS CAPSULORHEXIS FORCEPS ROUND BODIED STRAIGHT TIPPED |
| SC34M | CAPSULORHEXIS ROUND BODIED CURVED Pince Capsulorhexis Manche Rond / Courbe Kapsulorhexispinzette gebogen runder Griff |
| SC29M | MICS CAPSULORHEXIS FORCEPS CURVED Pince Capsulorhexis MICS Courbe / MICS Kapsulorhexispinzette gebogen |

FIXATION FORCEPS / PINCES / PINZETTEN - Pg10-12

| | |
|--------------|---|
| SC59 | ADSONS DISSECTION FORCEPS / Pince d'Adson à Griffes / Adson Dissektionspinzette |
| SC56 | BISHOPS FORCEPS / Pince de Bishop / Bishops Pinzette |
| SC38 | BONN FORCEPS / Pince de Bonn/ à Griffes/ Manche court / Irispinzette/ Bonn/ kurzer Griff |
| SC129 | BONN FORCEPS / Pince de Bonn/ à Griffes/ Manche long / Irispinzette/ Bonn/ langer Griff |
| SC22 | COLIBRI TOOTHED FORCEPS / Pince Colibri/ à Griffes/ Manche long / Kolibri Pinzette/ gezahnt/ langer Griff |
| SC20 | COLIBRI TYING FORCEPS / Pince Colibri/ à Griffes/ Manche court / Kolibri Pinzette/ gezahnt/ kurzer Griff |
| SC115 | FINE CLAW FORCEPS / Pince à Griffes Fines / feine Greifpinzette |
| SC55 | HOSKINS FORCEPS / Pince de Hoskins / Hoskins Pinzette |
| SC116 | IRIS FORCEPS / Pince à iris/ griffes fines / Irispinzette/ gezahnt/ sehr fein |
| SC128 | IRIS FORCEPS POINTED / Pince à iris/ Extrémités pointues/ Mors striés / Irispinzette/ gerade/ gerieft |
| SC117 | JAYLES FORCEPS / Pince de Jayles / Jayles Pinzette |
| SC49 | STRAIGHT NOTCHED FORCEPS / Pince à fixer alvéolée de Hoskins / Hoskins Pinzette gerade |
| SC54 | ST. MARTINS FORCEPS / Pince de Paufigue / St. Martins Pinzette |
| SC21 | STRAIGHT TOOTHED FORCEPS FINE / Pince à Tissus, droite / Pinzette gerade, fein mit Zähnen |

GENERAL FORCEPS / PINCES / PINZETTEN - Pg13-15

| | |
|--------------|---|
| SC119 | ADSONS FORCEPS / Pince d'Adson, Mors striés / Adson anatomische Pinzette, gerieft |
| SC113 | ADSONS FORCEPS FINE SERRATED / Pince d'Adson, Striée, Fine / Adson anatomische Pinzette, gerieft, sehr fein |

| | |
|--------------|--|
| SC104 | DENIS SCLERECTOMY FORCEPS / Pince à sclérectomie de Denis / Denis Sclerectomy Pinzette |
| SC58 | EPILATION - CILIA FORCEPS / Pince à Cil / Fremdkörper & Zilienpinzette |
| SC127 | EPILATION FORCEPS / Pince à épiler / Zilien-und Fremdkörperpinzette |
| SC60 | JEWELLERS FORCEPS / Pince de Joailler / Feine Fadenpinzette (Jewellers) |
| SA12 | MOORFIELDS FORCEPS / Pince de Moorfield / Moorfield Pinzette |
| SA17 | MOORFIELDS FORCEPS FINE / Pince de Moorfield Pointes fines / Moorfield Pinzette,fein |
| SG02 | MOSQUITO FORCEPS / Pince Mosquito / Mosquito Pinzette |

LENS FORCEPS / PINCES / PINZETTEN - Pg16

| | |
|--------------|---|
| SC112 | IOL REMOVAL FORCEPS / Pince à extraction d'implant / IOL Explantations Pinzette |
| SC112 | IOL REMOVAL FORCEPS |
| SC51 | LENS FOLDING FORCEPS / Pince à implants / IOL Faltpinzette |
| SC52 | LENS IMPLANTATION FORCEPS / Pince à implantation IOL / IOL Implantations Pinzette |

TYING FORCEPS / PINCES / PINZETTEN - Pg16

| | |
|--------------|---|
| SC50 | CURVED NOTCHED FORCEPS / Pince de Troutman / gebogene chirurgische Pinzette |
| SC73 | KELMAN MC PHERSON FORCEPS / Pince Kelman McPherson / Kelman Mc Pherson Pinzette |
| SC72 | MC PHERSON TYING FORCEPS / Pince McPherson / Angewinkelte Fadenpinzette (McPherson) |
| SC19 | STRAIGHT TYING FORCEPS / Pince à fil droite, MC / Knüpf-/Fadenpinzette gerade |
| SC71 | STRAIGHT TYING FORCEPS / Pince à fil droite, ML / Knüpf-/Fadenpinzette gerade |
| SC120 | SUTURE REMOVAL FORCEPS CURVED / Pince à sutures, courbe / Fadenpinzette zur Nahtentfernung, gebogen |
| SC105 | TROUTMAN TYING FORCEPS CURVED / Pince monofilament de Troutman / Knüpfpinzette nach Troutman, gebogen |

MANIPULATORS & CHOPPERS / MANIPULATEURS & CHOPPERS / MANIPULATEUREN & CHOPPER Pg19-23

| | |
|--------------|--|
| SC39 | CHOPPER AND DIALLER / Double Manipulateur Chopper/Y / Chopper und Spalter |
| SC95 | CHOPPER AND KLUGEN DOUBLE ENDED |
| SC97 | DENIS CLEAVER / Cliver de Denis / Denis Spalter |
| SC89 | DRYSDALE MICRO MANIPULATOR |
| SC47 | DRYSDALE NUCLEUS SPATULA |
| SC100 | GOVAN MANIPULATOR / Manipulateur Govan / Govan Manipulator |
| SC102 | GOVAN MANIPULATOR AND PHACO CHOPPER / Chopper Phaco & Manipulateur Govan / Govan Manipulator & Phako Chopper |
| SC45 | GREENS CHOPPER / Chopper de Greens / Greens Chopper |
| SC63 | IRIS REPOSITOR / Manipulateur à iris / Iris Repositor |
| SC90 | JONES CHOPPER / Chopper Jones / Jones Chopper |
| SC44 | KOCH CHOPPER / Chopper Koch / Koch Chopper |
| SC143 | KOCH PHACO SPATULA / Spatule de Koch / Phako Spatel nach Koch |
| SC106 | LE BUISSON MANIPULATOR / Manipulateur Le Buisson / Le Buisson Manipulator |
| SC46 | LENS DIALER (SINSKEY HOOK) / Crochet de Sinskey / Linsenhäkchen nach Sinskey |

- SC146** LESTER HOOK / Crochet de Lester / Lester Haken
SC48 MUSHROOM NUCLEUS MANIPULATOR / Manipulateur champignon / Mushroom Nukleus Manipulator
SC147 PHACO CHOPPER AND MUSHROOM MANIPULATOR DOUBLE ENDED / Manipulateur Double Chopper-
 Champignon / Phaco Chopper und Mushroom Manipulator (Kombi-Instrument)
SC43 PHACO CHOPPER / Chopper Phaco / Phako Chopper
SC42 PHACO CHOPPER CLAW / Chopper Claw, fin / Phako Claw Chopper (krallen förmig, fein)
SC101 PISACANO ROTATOR AND NEUHANN MANIPULATOR DOUBLE END / Manipulateur double Pisacano-
 Neuhann / Pisacano-Neuhann Manipulator (Kombi-instrument)
SC144 SINSKEY HOOK AND GOVAN MANIPULATOR DOUBLE ENDED / Manipulateur double Sinskey-Govan /
 Sinskey Häkchen und Govan Manipulator (Kombi-Instrument)
SC41 Y SHAPED NUCLEUS ROTATOR / Manipulateur de Noyau Y, Vertical / Y-Nukleus Rotator

SPECULUMS / BLEPHAROSTATS / LIDSPERRER - Pg25-27

- SC01** BROAD GAUGE WIRE SPECULUM / Blépharostat fils large / Breiter Gauge Draht Lidsperrerr
SC124 HEAVY GAUGE WIRE SPECULUM OPEN LOOP / Blépharostat fil, grand diamètre / Drahtlidsperrerr groß,
 offen
SC09 KRATZ BARRAQUER CLOSED LOOP / Blépharostat de Barraquer, fenêtré / Lidsperrerr nach Kratz,
 gefenstert
SC132 KRATZ BARRAQUER CLOSED LOOP SPECULUM / Blépharostat Kratz-Barraquer, fenêtré / Lidsperrerr Kratz-
 Barraquer, gefenstert
SC11 KRATZ BARRAQUER OPEN LOOP WIRE SPECULUM
SC107 KRATZ OPEN LOOP WIRE SPECULUM (BLUE) / Blepharostat de Kratz (bleu) / Kratz Lidsperrerr mit offenen
 Valven (Blau)
SC05 KRATZ SOLID BLADED WIRE SPECULUM / Blépharostat Kratz-Barraquer, anses pleines (Chrome) / Kratz-
 Barraquer Lidsperrerr, Valven aus Chrom (solid)
SC13 KRATZ CLOSED LOOP SPECULUM / Blépharostat temporal de Barraquer / Lidsperrerr nach Barraquer,
 feiner gebogener Draht temporal offen
SC150 OPEN LOOP WIRE SPECULUM
SC155 SPECULUM OPEN LOOP
 TEMPORAL SPECULUM / Blépharostat temporal de Kratz / Lidsperrerr nach Kratz, feiner gebogener
SC12 Draht temporal offen
 WIDE ANGLED CLOSED LOOP SPECULUM
SC07 CLOSED LOOP WIDE SPECULUM / Blepharostat anses fils large / Spekulum/Lidsperrerr mit
SC07 geschlossenen Valven

SPECULUMS MECHANICAL - Pg27-28

- SC03** BANGERTER ADULT SPECULUM / Blepharostat de Bangertter / Lidsperrerr für Erwachsene nach Bangertter
SR10 BLADED ADJUSTABLE SPECULUM / Blépharostat mécanique, anses pleines / Lidsperrerr nach Barraquer,
 Valven aus Vollmaterial
SR11 KRATZ SCREW ADJUSTABLE WIRE SPECULUM OPEN LOOP / Blépharostat mécanique Kratz, anses fil /
 mechanisches verstellbarer Lidsperrerr nach Kratz (feiner gebogener Draht) temporal offen

- SC126** CLARKE SPECULUM / Blépharostat de Clarke / Clarke Lidsperrerr
SC121 Lieberman Speculum
SC08 SOLID BLADED SCREW ADJUSTABLE SPECULUM / Blépharostat mécanique / Mechanisch feststellbarer
 Lidsperrerr mit Valven aus Vollmaterial

PAEDIATRIC SPECULUM - Pg29

- SC10** KRATZ PAEDIATRIC CLOSED LOOP 10 MM SPECULUM / Blépharostat pédiatrique de Kratz /
 paediatrischer Lidsperrerr nach Kratz 10mm
SC103 NEONATE PAEDIATRIC SPECULUM / Blépharostat néonatal Necker / Neonatal paediatrischer Lidsperrerr
SC15 NEONATE 4 MM SPECULUM / Blépharostat néonatal 4mm / Neonatal Lidsperrerr 4mm
SC14 PAEDIATRIC 6 MM SPECULUM / Blépharostat pédiatrique 6mm / paediatrischer Lidsperrerr 6mm

SCISSORS / CISEAUX / SCHEREN - Pg31-33

- SC82/83** CASTROVIEJO CORNEAL SCISSORS RIGHT AND LEFT / Ciseaux de Castroviejo droite et gauche /
 Castroviejo Cornea Schere rechts und links
SC111 CASTROVIEJO SCISSORS STRAIGHT / Ciseaux de Castroviejo, droits / Castroviejo Schere, gerade
SC110 CASTROVIEJO SCISSORS STRAIGHT & ROUND / Ciseaux de Castroviejo, courbes, mousses / Castroviejo
 Mikro-Schere, gebogen, stumpf
SC80 ONG IRIS SCISSORS / Ciseaux ONG / ONG Iris Schere
SC109 OSHER SCISSORS / Ciseaux à implant (Osher) / Osher Schere
SA13 STALLARD CORNEAL SCISSORS / Ciseaux de Stallard / Stallard Hornhautschere
SC81 VANNAS SCISSORS / Ciseaux de Vannas / Vannas Schere
SC108 VANNAS SCISSORS CURVED / Ciseaux de Vannas, courbes / Vannas Schere gebogen
SA11 WESTCOTT CONJUNCTIVAL SCISSORS / ciseaux de Westcott, courbes / Westcott Bindeha
 u t s c h e r e
SA14 WESTCOTT SCISSORS DOLPHIN NOSE / Ciseaux de Westcott, extrémités émoussées (en forme de
 dauphin) / Westcott Schere, stumpfe Enden (delphinförmig)
SA15 WESTCOTT DOLPHIN NOSE SCISSORS BLUNT TIPPED

SCISSORS GENERAL / CISEAUX / SCHEREN - Pg33-35

- SG35** DRESSING SCISSORS / Ciseaux à disséquer, extrémités courbes et mousses / Verbandschere gebogen
 mit abgerundete Spitze
SG34 ENUCLEATION SCISSORS / Ciseaux à énucléation / Enucléationsschere
SG14 GENERAL LONG HANDLE SCISSORS BLUNT END / Ciseaux courbes, bout mousse, manche long /
 Augenschere lang, stumpfe Enden
SG06 GENERAL SCISSORS BLUNT END / Ciseaux courts, pointes droites, bout mousse / Augenschere, stumpfe
 Enden
SG01 GENERAL SCISSORS SHARP POINTED / Ciseaux courts pointus / Augenschere, spitze Enden
SC96 IRIS SCISSOR
SG30 METZENBAUM SCISSORS CURVED POINTED / Ciseaux de Metzenbaum, courbes, pointus / Metzebaum
 Schere, gebogen, spitz

- SG07** STRABISMUS SCISSORS / Ciseaux à strabisme / Strabismus Schere
SG23 STEVENS SCISSORS CURVED / Ciseaux de Stevens courbes / Schere nach Stevens, gebogen
SG31 STEVENS TENOTOMY SCISSORS POINTED

GENERAL OPHTHALMIC / OPHTHALMOLOGIE GENERALE/ ALLGEMEINE OPHTHALMOLOGIE - Pg37-40

- SG32** BALL TIP MUSCLE HOOK / Crochet à muscle, point boule / Schielhäkchen, ballförmige Spitze
SG92 BULLDOG CLIP / Bulldog Clip / Bulldog Clip
SC91 CHAVASSE HOOK / Crochet de Chavasse / Chavasse Häkchen
SG20 DESMARRES RETRACTOR / Retracteur de Desmarre / Desmarres Retraktor
SG17 EYELID PLATE / Plaque à paupière de Jaeger / Lidplatte
SG28 GILLIES HOOK SMALL
SG37 GUTHRIE HOOK / Crochet de Guthrie / Guthrie Häkchen
- SG33** HALSTEAD MOSQUITO FORCEPS STRAIGHT / Pince de Halstead, droite / Halstead Arterienklemme, gerade
SG10 JAMESON STRABISMUS HOOK 5 MM / Crochet à strabisme de Jameson / Jameson Schielhaken
SG15 KUGLEN HOOK
SG22 MOSQUITO FORCEPS CURVED / Pince Mosquito, courbe / Mosquito Pinzette, gebogen
SC23 SCLERAL DEPRESSOR BULBOUS END / Dépresseur scléral / Skleral Depressor wulstige Enden
SC25 SCLERAL DEPRESSOR DOUBLE ENDED / Dépresseur scléral double de Schocket / Skleral Depressor zweiseitig
SC24 SCLERAL DEPRESSOR HAMMER END / Dépresseur scléral courbe / Skleral Depressor hammerförmiges Ende
SG12 THORNTON FIXATION RING
SG13 THORNTON FIXATION RING STANDARD 16 MM / Anneau de fixation de Thornton 16mm / Thornton Fixationsring standard 16mm
SG03 TRABECULECTOMY PUNCH / Punch à trabéculectomie / Trabekulotomie

CORNEAL & REFRACTIVE / CORNEE & REFRACTIVE / KORNEAL & REFRAKTIV - Pg43-50

- SR41** ALCOHOL WELL VARIOUS SIZES / Anneau marqueur / Lasek Alkohol Trichter
SR42 ALCOHOL WELL VARIOUS SIZES / Anneau marqueur / Lasek Alkohol Trichter
SR43 ALCOHOL WELL VARIOUS SIZES / Anneau marqueur / Lasek Alkohol Trichter
SR109 DMEK BERNHEIM HOOK / Crochet DMEK de Bernheim / DMEK Häkchen nach Bernheim
SR100 CORNEAL EPITHELIUM SPATULA ROUND / Spatule pour épithélium cornéen, ronde / Epithélium Hornhautspatel, rund
SR72 CORNEAL MARKER 8 LINE
SR40 CK LASIK MARKER / Marqueur LASIK avec croix de centrage / Hornhautmarker mit Fadenkreuz
SR04 CORNEAL RING PUSHER / Manipulateur anneau intra-cornéen / Manipulator für Cornealen Ring
SR05 CORNEAL RING POCKETING HOOK / Crochet anneau intra-cornéen / Taschen Haken für Cornealen Ring

- SR06** CORNEAL RING POCKET SPREADER / Spatule anneau intra-cornéen / Taschen Spreizer für Cornealen Ring
SR07 CORNEAL RING FORCEPS KERARING / Pince à anneau intra-cornéen - Keraring / Keraring Pinzette für Cornealen Ring
SR08 CORNEAL RING FORCEPS INTACS / Pince à anneau intra-cornéen - INTACS / Intacs Pinzette für Cornealen Ring
SR61 EPI PEELER RECTANGLE SHAPED / Scarificateur, forme D / Epi Peeler - quadratische Form (LASIK)
SR60 EPI PEELER D SHAPED / Scarificateur Micro-Houe / Epi Peeler, gebogen
SR51 EPI PEELER / Scarificateur cornéen / Epi Peeler
SR91 IMPLANTATION FORCEPS
SR34 INTRA LASIK FLAP SEPARATOR / Elévateur volet cornéen, crochet pointu / Intra-Lasik Flap Separator
SR102 KOCABA DMEK FORCEPS / P i n c e D M E K d e K o c a b a / K o c a b a D M E K P i n z e t t e
SR21 LASIK CORNEAL MARKER / Marqueur cornéen en fourche / Lasik Hornhautmarker
SR02 LASIK FEMTO SEPERATOR AND ELEVATOR DOUBLE ENDED / Séparateur-Élévateur Siebel / Lasik Femto Seperator und Elevator zweiseitig
SR31 LASIK FLAP ELEVATOR / Spatule Flap Lasik / LASIK Flap Heber
SR14 LASIK FLAP-LIFTING FORCEPS / Pince angulée mors long / Lasik Flap Hebebinzette
SR71 LASIK FLAP FORCEPS / Pince Flap Lasik / Lasik Flap Pinzette
SR103 LASIK FLAP SEPERATOR AND ELEVATOR DOUBLE END / Séparateur-élévateur LASIK, droit / Lasik Flap Separator / Flap Heber gerade
SR01 LENTICULE EXTRACTION FORECEPS
SR108 PATON SPATULA & SPOON DOUBLE ENDED / Spatule double de Paton / Paton Spatel und Löffel Zweiseitig
SR104 REINSTEIN LENTICULE SEPERATOR HOOK AND SPATULA DOUBLE END / Relex Crochet / Séparateur, Spatule / Relex Häkchen / Separator (Kombi-Instrument)
SR70 REINSTEIN LENTICULE SEPERATOR
SR101 S MARKER / Marqueur DMEK «S» / S Marker
SR03 SINSKEY AND REVERSE SINSKEY HOOK / Manipulateur double Sinskey-Sinskey inversé / Sinskey Häkchen umgekehrt Doppelinstrument
SR105 TREPHINE VARIOUS SIZES / Trépan cornéen / Trepane verschiedene Größen
SR106 TREPHINE VARIOUS SIZES / Trépan cornéen / Trepane verschiedene Größen
SR107 TREPHINE VARIOUS SIZES / Trépan cornéen / Trepane verschiedene Größen
SR110 TREPHINE VARIOUS SIZES / Trépan cornéen / Trepane verschiedene Größen

MEASURERS & MARKERS / MESUREURS & MARQUEURS / MESSINSTRUMENTE & MARKER - Pg53-54

- SG04** BRAUNSTEIN FIXED CALIPER 1.8 MM THICKNESS / Marqueur de Braunstein, métal / Braunstein Kaliber fix 1.8mm Dicke
SG05 CASTROVIEJO ADJUSTABLE CALIPER / Compas de Castroviejo, pointes droites / Castroviejo verstellbarer Zirkel
SC130 CORNEAL MARKER GUIDE
SG93 HELVESTON RULER / Règle d'Helveston / Helveston Meßspatel

- SC65** MENDEZ GAUGE / Anneau de Mendez / Messlehre nach Mendez (zur Bestimmung des Einschnittwinkels)
- SG19** SCLERAL MARKER ROUND BODIED DOUBLE ENDED / Marqueur scléral métal 3.5 mm - 4.5 mm / Skleral Marker runder Griff, zweiseitig

CURETTES & CLAMPS / CURETTES & PINCES / KURETTEN & KLEMMEN - Pg57

- CR01**
- CR02**
- CR03** CURETTE VARIOUS SIZES / Curettes à chalazion / Küretten verschieden Größen
- CR05** CURETTE DOUBLE ENDED / Curette double / Kürette zweiseitig
- CS01** CHALAZION CLAMP SMALL / Pince à chalazion 12mm / Chalazion Klemme klein
- CS02** CHALAZION CLAMP LARGE / Pince à chalazion large, ovale / Chalazion Klemme groß

NEEDLE HOLDERS / PORTE-AIGUILLES / NADELHALTER - Pg59

- SC142** BARRAQUER NEEDLE HOLDER CURVED / Porte-Aiguille de Barraquer, courbe / Barraquer Nadelhalter gebogen
- SC118** CASTROVIEJO NEEDLE HOLDER LOCKING / Porte-aiguille de Castroviejo droit, verrou / Castroviejo Nadelhalter mit Sperre
- SC62** NEEDLE HOLDER BLUNT LOCK / Porte-aiguille, bout mousse, verrou / Nadelhalter stumpf, mit Sperre
- SC114** NEEDLE HOLDER CURVED / Porte-aiguille, courbe / Nadelhalter, gebogen
- SC61** SUTURE NEEDLE HOLDER CURVED

LACRIMAL / LACRYMAL / LAKRIMAL - Pg61-62

- PD05** BOWMANS PROBE SIZE 00-0 DOUBLE ENDED / Sonde de Bowman double 00-0 / Tränenwegsonde Bowman Größe 00-0 zweiseitig
- PD06** BOWMANS PROBE SIZE 000-00 DOUBLE ENDED / Sonde de Bowman double 000-00 / Tränenwegsonde Bowman Größe 000-00 zweiseitig
- PD08** LACRIMAL PROBE / Sonde à voies lacrymales / Lakrimal Sonde
- PD02** PUNCTAL DILATOR BLUNT / Dilatateur voies lacrymales, mousse / Tränenwegsonde stumpf
- PD03** PUNCTAL DILATOR DOUBLE ENDED / Dilatateur voies lacrymales, fin-medium / Tränenwegsonde zweiseitig
- PD01** PUNCTAL DILATOR SHARP / Dilatateur voies lacrymales, très fin / Tränenwegsonde scharf
- PD04** PUNCTAL DILATOR SHORT / Dilatateur voies lacrymales, manche court / Tränenwegsonde kurz
- PD07** PUNCTAL DILATOR & LACRIMAL PROBE DOUBLE ENDED / Dilatateur-Sonde pour voies lacrymales / Tränenwegsonde zweiseitig





IQ MEDICAL – Exclusive Australian distributors

+61 8 8357 8022 | sales@iqmedical.com.au

2/86 Mary Street, Unley SA 5061 | iqmedical.com.au



M-0050-A