

# TrabEx PRO™ – goniotomy, refined



## Confidence through visibility, with comfort and control

MicroSurgical Technology's TrabEx PRO™ is an innovative dual-blade excisional goniotomy device designed to support controlled partial excision of the trabecular meshwork during MIGS procedures. The uniquely integrated active irrigation and aspiration enable stable chamber control and clear visualisation throughout the procedure.

### FEATURES AND BENEFITS

- Dual serrated trapezoidal blades for controlled excision of trabecular meshwork.
- Active irrigation and aspiration provide a differentiated level of chamber stability and clear visualisation.
- Rounded footplate designed to accommodate varying patient anatomies.
- Flexible silicone sleeve enabling compatibility with 2.2–2.4 mm incisions.
- Sterile, single use design for procedural efficiency and consistency.

Access the latest product  
information here



**IQ MEDICAL – Exclusive Australian distributors**

+61 8 8357 8022 | [sales@iqmedical.com.au](mailto:sales@iqmedical.com.au)

2/86 Mary Street, Unley SA 5061 | [iqmedical.com.au](http://iqmedical.com.au)

M-0090-A

