

Packer/Chang IOL cutting system

MST Complex Cases,
Simplified™



A single-use system exclusive to IQ Medical

MicroSurgical Technology's **Packer/Chang IOL cutting system** is the gold standard for safe, controlled IOL explantation. The combination box includes 3 single-use Packer/Chang IOL cutters and 3 single-use micro-holding forceps. Designed for firm grasping and precise cutting through small incisions, it simplifies complex cases with confidence.

FEATURES AND BENEFITS

- Low-profile 19G Packer/Chang IOL cutters for minimal corneal distortion.
- Hingeless design reduces lens tilt during cutting.
- Cuts all foldable IOL materials, including acrylic and silicone.
- 23G serrated micro-holding forceps for secure IOL handling.
- Supplied sterile in a combination box for convenience.



Access the MST
catalogue here



IQ MEDICAL – Exclusive Australian distributors

+61 8 8357 8022 | sales@iqmedical.com.au

2/86 Mary Street, Unley SA 5061 | iqmedical.com.au

M-0108-A

